

Anti-Uqcc2 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9BrT2

Item Number: CM0027353

Introduction

Rabbit polyclonal antibody to human UQCC2, validated for WB, IHC-P, ELISA. Detects human, mouse, rat UQCC2; involved in mitochondrial complex III assembly and respiratory chain function.

[Learn More](#)

Parameter	Value
Product Name	UQCC2 Rabbit pAb
Remarks/Alias	M19; Cbp6; MNF1; MC3DN7; C6orf125; C6orf126; bA6B20.2; UQCC2
Species	Human
Gene ID (Human)	84300
Gene ID	84300
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 14-126 of human UQCC2 (NP_115716.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BRT2
Protein Weight	15kDa
Shipping	Ice bag