


Anti-Uqcr10 Polyclonal Antibody For Wb And Elisa - Q9Udw1

Item Number: CM0029781



Introduction

Rabbit polyclonal antibody targeting human UQCR10 (Cytochrome b-c1 complex subunit 9), ideal for Western blot and ELISA. Detects human, mouse, and rat samples. ~7 kDa mitochondrial inner membrane protein of Complex III.

[Learn More](#)

Parameter	Value
Product Name	UQCR10 Rabbit pAb
Remarks/Alias	QCR9; UCRC; HSPC051; HSPC119; HSPC151; UCCR7.2; UQCR10
Species	Human
GeneID (Human)	29796
GeneID	29796
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-63 of human UQCR10 (NP_037519.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UDW1
Protein Weight	7kDa
Shipping	Ice bag