


## Anti-Uqcr10 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Udw1

Item Number: CM0016399



### Introduction

Rabbit polyclonal antibody targeting human UQCR10 (Complex III subunit 9). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for mitochondrial electron transport chain research.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | UQCR10 Rabbit pAb                                                                                                |
| Remarks/Alias    | QCR9, UCRC, HSPC051, HSPC119, HSPC151, UCCR7.2, UQCR10                                                           |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 29796                                                                                                            |
| GeneID           | 29796                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-63 of human UQCR10 (NP_037519.2) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | Q9UDW1                                                                                                           |
| Protein Weight   | 7kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                          |