


Anti-Uqcrb Rabbit Polyclonal Antibody For Wb, Elisa - P14927

Item Number: CM0017640



Introduction

UQCRB Rabbit Polyclonal Antibody validated for WB, ELISA in human, mouse, rat. Targets mitochondrial complex III subunit 7, driving electron transport and oxidative phosphorylation. Supports mitochondrial disease and respiratory chain research.

[Learn More](#)

Parameter	Value
Product Name	UQCRB Rabbit pAb
Remarks/Alias	QPC; QCR7; QP-C; UQBC; UQBP; UQPC; UQCR6; MC3DN3; UQCRB
Species	Human
GeneID (Human)	7381
GeneID	7381
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-111 of human UQCRB (NP_006285.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14927
Protein Weight	14kDa
Shipping	Ice bag