


Anti-Uqcrc1 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P31930

Item Number: CM0004850



Introduction

High-performance rabbit monoclonal antibody against UQCRC1 (Complex III subunit 1), validated for WB, IHC-P, IF/ICC, IP, and ELISA. Suitable for human, mouse, and rat samples. Ideal for mitochondrial respiration and Parkinson's disease research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] UQCRC1 Rabbit mAb
Remarks/Alias	QCR1; PKNPY; UQCR1; D3S3191
Species	Human
Gene ID (Human)	7384
Gene ID	7384
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31930
Protein Weight	53kDa
Shipping	Ice bag