


Anti-Uqcrc1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P31930

Item Number: CM0026513



Introduction

High-quality rabbit polyclonal antibody against human UQCRC1, validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for mitochondrial research and complex III studies.

[Learn More](#)

Parameter	Value
Product Name	UQCRC1 Rabbit pAb
Remarks/Alias	QCR1; PKNPY; UQCR1; D3S3191; UQCRC1
Species	Human
GeneID (human)	7384
GeneID	7384
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 156-480 of human UQCRC1 (NP_003356.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31930
Protein Weight	53kDa
Shipping	Ice bag