


Anti-Uqcrc2 Polyclonal Antibody (Rabbit) For Wb, Ihc-P, If/Icc, Ip, Elisa - P22695

Item Number: CM0030736



Introduction

High-quality rabbit polyclonal antibody targeting human UQCRC2, validated for WB, IHC-P, IF/ICC, IP, and ELISA. Cross-reacts with mouse and rat. Ideal for mitochondrial complex III and electron transport chain research.

[Learn More](#)

Parameter	Value
Product Name	UQCRC2 Rabbit pAb
Remarks/Alias	QCR2; UQCR2; MC3DN5; UQCRC2
Species	Human
GeneID (Human)	7385
GeneID	7385
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 224-453 of human UQCRC2 (NP_003357.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22695
Protein Weight	48kDa
Shipping	Ice bag