


Anti-Uqcrc2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P22695

Item Number: CM0011395



Introduction

UQCRC2 Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA. Targets human UQCRC2 (P22695), a mitochondrial complex III subunit essential for electron transport and ATP synthesis. Ideal for mitochondrial disease research.

[Learn More](#)

Parameter	Value
Product Name	UQCRC2 Rabbit mAb
Remarks/Alias	QCR2; UQCR2; MC3DN5; UQCRC2
Species	Human
GeneID (Human)	7385
GeneID	7385
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22695
Protein Weight	48kDa
Shipping	Ice bag