


Anti-Uqcrfs1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P47985

Item Number: CM0017877



Introduction

Rabbit polyclonal antibody against human UQCRFS1 (Rieske iron-sulfur protein), validated for WB, IHC-P, ELISA. Detects endogenous UQCRFS1 and cross-reacts with mouse. Suitable for mitochondrial complex III research.

[Learn More](#)

Parameter	Value
Product Name	UQCRFS1 Rabbit pAb
Remarks/Alias	RIP1; RIS1; RISP; UQCR5; MC3DN10; UQCRFS1
Species	Human
Gene ID (Human)	7386
Gene ID	7386
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 79-274 of human UQCRFS1 (NP_005994.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	P47985
Protein Weight	30kDa
Shipping	Ice bag