


Anti-Uqcrfs1 (Rieske Iron-Sulfur Protein) Polyclonal Antibody For Wb, Ihc-P, Elisa - P47985

Item Number: CM0030976



Introduction

Rabbit polyclonal antibody against human UQCRFS1 (Rieske iron-sulfur protein), a mitochondrial Complex III subunit. Validated for Western blot, IHC-P, and ELISA; reactive with human, mouse, and rat samples. Ideal for oxidative phosphorylation research.

[Learn More](#)

Parameter	Value
Product Name	UQCRFS1 Rabbit pAb
Remarks/Alias	RIP1, RIS1, RISP, UQCR5, MC3DN10, UQCRFS1
Species	Human
GeneID (Human)	7386
GeneID	7386
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 141-274 of human UQCRFS1 (NP_005994.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P47985
Protein Weight	30kDa
Shipping	Ice bag