


Anti-Uqcrcq Polyclonal Antibody For Wb, Elisa - O14949

Item Number: CM0025111



Introduction

High-quality rabbit anti-UQCRCQ polyclonal antibody for Western blot and ELISA, recognizing human and mouse samples. Targets complex III subunit 8 essential for mitochondrial electron transport and oxidative phosphorylation.

[Learn More](#)

Parameter	Value
Product Name	UQCRCQ Rabbit pAb
Remarks/Alias	QPC; QCR8; QP-C; UQCR7; MC3DN4; UQCRCQ
Species	Human
Gene ID	27089
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-82 of human UQCRCQ (NP_055217.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O14949
Protein Weight	10kDa
Shipping	Ice bag