


Anti-Uridine-Cytidine Kinase 2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bzx2

Item Number: CM0013392



Introduction

Rabbit polyclonal antibody targeting human UCK2, validated for WB and ELISA. Cross-reacts with mouse UCK2. UCK2 is a uridine-cytidine kinase involved in nucleoside salvage pathways.

[Learn More](#)

Parameter	Value
Product Name	UCK2 Rabbit pAb
Remarks/Alias	UK; UMPK; TSA903; UCK2
Species	Human
Gene ID (Human)	7371
Gene ID	7371
Immunogen	Recombinant protein - Recombinant fusion protein containing a sequence corresponding to amino acids 1-261 of human UCK2 (NP_036606.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BZX2
Protein Weight	29kDa
Shipping	Ice bag