


Anti-Urocortin-2 (Ucn2) Rabbit Polyclonal Antibody For Wb And Elisa - Q96Rp3

Item Number: CM0029950



Introduction

High-quality rabbit polyclonal antibody targeting human UCN2, validated for Western blot and ELISA. Detects human, mouse, and rat. Predicted MW 12 kDa. Ideal for stress response, cardiovascular, and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	UCN2 Rabbit pAb
Remarks/Alias	UR; SRP; URP; UCN1; UCN-II; UCN2
Species	Human
Gene ID (Human)	90226
Gene ID	90226
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-112 of human UCN2 (NP_149976.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96RP3
Protein Weight	12kDa
Shipping	Ice bag