

Anti-Urod Polyclonal Antibody For Wb, Elisa - P06132

Item Number: CM0025455



Introduction

Rabbit polyclonal antibody against human UROD (uroporphyrinogen decarboxylase), applicable in WB and ELISA. Cross-reacts with mouse. Ideal for heme biosynthesis and porphyria research. 41 kDa.

[Learn More](#)

Parameter	Value
Product Name	UROD Rabbit pAb
Remarks/Alias	PCT, UPD, UROD
Species	Human
GeneID	7389
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-367 of human UROD (NP_000365.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P06132
Protein Weight	41kDa
Shipping	Ice bag