

Anti-Uroplakin Iii/Upk3A Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q9Jkx8

Item Number: CM0028591

Introduction

Rabbit polyclonal antibody against Uroplakin III (UPK3A), suited for WB, IF/ICC, IHC-P, and ELISA. Detects mouse, human, and rat UPK3A. Ideal for urothelial biology and bladder disease research.

[Learn More](#)

Parameter	Value
Product Name	Uroplakin III/UPK3A Rabbit pAb
Remarks/Alias	UP3A; UPK3; UPIII; UPIIIA
Species	Mouse
Gene ID (Human)	7380
Gene ID	22270
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 188-287 of mouse Uroplakin III/UPK3A (NP_075967.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9JKX8
Protein Weight	31 kDa
Shipping	Ice bag