

Anti-Usf1 Rabbit Monoclonal Antibody For Wb And Elisa - P22415

Item Number: CM0002018



Introduction

Rabbit monoclonal antibody targeting USF1 (Upstream stimulatory factor 1), suitable for WB and ELISA applications. Reacts with human and mouse USF1. Ideal for transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	USF1 Rabbit mAb
Remarks/Alias	UEF; FCHL; MLTF; FCHL1; MLTF1; HYPLIP1; bHLHb11; USF1
Species	Human
GeneID (human)	7391
GeneID	7391
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human USF1 (P22415)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P22415
Protein Weight	34kDa
Shipping	Ice bag