

Anti-Usf1 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P22415

Item Number: CM0006957



Introduction

High-quality rabbit monoclonal antibody targeting human USF1 (SWISS P22415). Validated for WB, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat. Suitable for transcription research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] USF1 Rabbit mAb
Remarks/Alias	UEF; FCHL; MLTF; FCHL1; MLTF1; HYPLIP1; bHLHb11; [KD Validated] USF1
Species	Human
GeneID (Human)	7391
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-310 of human USF1. (NP_009053.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22415
Protein Weight	34kDa
Shipping	Ice bag