

Anti-Usf1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P22415

Item Number: CM0016502



Introduction

USF1 rabbit polyclonal antibody for Western blot, IHC-P, and ELISA applications. Detects USF1 in human, mouse, and rat samples. Targets the E-box-binding transcription factor involved in gene regulation.

[Learn More](#)

Parameter	Value
Product Name	USF1 Rabbit pAb
Remarks/Alias	UEF; FCHL; MLTF; FCHL1; MLTF1; HYPLIP1; bHLHb11; USF1
Species	Human
GeneID (human)	7391
GeneID	7391
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-310 of human USF1 (NP_009053.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22415
Protein Weight	34kDa
Shipping	Ice bag