


Anti-Usf2 Rabbit Polyclonal Antibody For Wb, Elisa - Q15853

Item Number: CM0026647



Introduction

Rabbit polyclonal antibody against human USF2 (FIP/bHLHb12) validated for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for transcription factor research and E-box binding studies.

[Learn More](#)

Parameter	Value
Product Name	USF2 Rabbit pAb
Remarks/Alias	FIP; bHLHb12; USF2
Species	Human
GeneID (human)	7392
GeneID	7392
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15853
Protein Weight	37kDa
Shipping	Ice bag