

Anti-Usol Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O60763

Item Number: CM0003079



Introduction

Rabbit monoclonal antibody against human USO1/p115, validated for WB, IHC-P, IF/ICC, and ELISA. Detects endogenous USO1 in human, mouse, and rat samples. Protein weight approximately 108 kDa. Ideal for Golgi transport research.

[Learn More](#)

Parameter	Value
Product Name	USO1 Rabbit mAb
Remarks/Alias	TAP; VDP; P115; USO1
Species	Human
Gene ID (Human)	8615
Gene ID	8615
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human USO1 (O60763).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60763
Protein Weight	108kDa
Shipping	Ice bag