

Anti-Usol Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O60763

Item Number: CM0021438



Introduction

Rabbit polyclonal antibody against USO1 (TAP/VDP/P115), validated for WB, IF/ICC, and ELISA. Reacts with human and rat. Ideal for Golgi vesicle transport studies.

[Learn More](#)

Parameter	Value
Product Name	USO1 Rabbit pAb
Remarks/Alias	TAP; VDP; P115; USO1
Species	Human
Gene ID (Human)	8615
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 663-962 of human USO1 (NP_003706.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	O60763
Protein Weight	108kDa
Shipping	Ice bag