

Anti-Usol Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O60763

Item Number: CM0027748



Introduction

Rabbit polyclonal antibody targeting human USO1 (p115/VDP/TAP), validated for WB, IHC-P, IF/ICC, and ELISA. Ideal for studying ER-Golgi vesicular transport and Golgi stack maintenance.

[Learn More](#)

Parameter	Value
Product Name	USO1 Rabbit pAb
Remarks/Alias	TAP; VDP; P115; USO1
Species	Human
GeneID (human)	8615
GeneID	8615
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1100-1200 of human USO1 (O60763).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O60763
Protein Weight	108kDa
Shipping	Ice bag