


Anti-Usp1 Rabbit Polyclonal Antibody For Wb - 094782

Item Number: CM0029219



Introduction

Rabbit polyclonal antibody targeting human USP1, validated for WB and ELISA. Detects USP1 in human and mouse samples. Useful for DNA damage repair and ubiquitin pathway research.

[Learn More](#)

Parameter	Value
Product Name	USP1 Rabbit pAb
Remarks/Alias	UBP; USP1
Species	Human
Gene ID (Human)	7398
Gene ID	7398
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human USP1 (NP_003359.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O94782
Protein Weight	88kDa
Shipping	Ice bag