

Anti-Usp11 Rabbit Monoclonal Antibody For Wb And Elisa - P51784

Item Number: CM0002491



Introduction

Rabbit monoclonal antibody against human USP11 (UHX1), validated for Western blotting (WB) and ELISA. Suitable for deubiquitination, NF-kappa-B signaling, DNA repair, and chromatin regulation studies.

[Learn More](#)

Parameter	Value
Product Name	USP11 Rabbit mAb
Remarks/Alias	UHX1; USP11
Species	Human
GeneID (human)	8237
GeneID	8237
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 800-900 of human USP11 (P51784).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51784
Protein Weight	110kDa
Shipping	Ice bag