

Anti-Usp11 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P51784

Item Number: CM0019781



Introduction

High-quality rabbit polyclonal antibody against USP11, validated for Western blot, IF/ICC, and ELISA. Ideal for ubiquitin pathway, DNA repair, and chromatin research across human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	USP11 Rabbit pAb
Remarks/Alias	UHX1; USP11
Species	Human
Gene ID	8237
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 70-190 of human USP11 (NP_004642.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51784
Protein Weight	110kDa
Shipping	Ice bag