

Anti-Usp13 Rabbit Polyclonal Antibody For Wb, Elisa - Q92995

Item Number: CM0014925



Introduction

Anti-USP13 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human, mouse, and rat USP13, a deubiquitinase regulating autophagy, ERAD, and cell cycle. Ideal for oncology and cell biology research.

[Learn More](#)

Parameter	Value
Product Name	USP13 Rabbit pAb
Remarks/Alias	ISOT3; IsoT-3; USP13
Species	Human
Gene ID (Human)	8975
Gene ID	8975
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 90-200 of human USP13 (NP_003931.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q92995
Protein Weight	97kDa
Shipping	Ice bag