


Anti-Usp14 Monoclonal Antibody For Wb, Ip, Elisa - P54578

Item Number: CM0003579



Introduction

Rabbit monoclonal antibody targeting USP14, a proteasome-associated deubiquitinase involved in ubiquitin recycling, CXCR4 degradation, and innate immunity. Validated for WB, IP, ELISA; cross-reacts with human, mouse, rat. Predicted MW 56kDa.

[Learn More](#)

Parameter	Value
Product Name	USP14 Rabbit mAb
Remarks/Alias	TGT; Ubp6; USP14
Species	Human
GeneID (Human)	9097
GeneID	9097
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 371-473 of human USP14 (P54578).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54578
Protein Weight	56kDa
Shipping	Ice bag