

Anti-Usp14 Polyclonal Antibody For Wb, Elisa - P54578

Item Number: CM0021287



Introduction

Rabbit polyclonal antibody targeting human USP14, validated for Western blot (WB) and ELISA. Suitable for research on ubiquitin-dependent proteolysis, innate immunity, and neuromuscular function.

[Learn More](#)

Parameter	Value
Product Name	USP14 Rabbit pAb
Remarks/Alias	TGT; Ubp6; USP14
Species	Human
GeneID (human)	9097
GeneID	9097
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 400-500 of human USP14 (NP_005142.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P54578
Protein Weight	56kDa
Shipping	Ice bag