

Anti-Usp15 Monoclonal Antibody For Wb, Ip, Elisa - Q9Y4E8

Item Number: CM0008420



Introduction

High-quality rabbit monoclonal antibody targeting USP15 (deubiquitinating enzyme 15). Ideal for studying TGF-beta, NF-kappa-B, and mitophagy pathways in human, mouse, rat tissues via WB, IP, ELISA.

[Learn More](#)

Parameter	Value
Product Name	USP15 Rabbit mAb
Remarks/Alias	UNPH4; UNPH-2
Species	Human
GeneID (Human)	9958
GeneID	9958
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 968-981 of human USP15 (NP_001239007.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4E8
Protein Weight	27 kDa/109-112 kDa
Shipping	Ice bag