

Anti-Usp15 Rabbit Polyclonal Antibody For Wb, Ip, If/Icc, Elisa - Q9Y4E8

Item Number: CM0013423



Introduction

Anti-USP15 Rabbit Polyclonal Antibody validated for WB, IP, IF/ICC, and ELISA. Detects human, mouse, rat USP15. Ideal for investigating deubiquitination, TGF-beta signaling, NF-kappa-B regulation, and mitophagy.

[Learn More](#)

Parameter	Value
Product Name	USP15 Rabbit pAb
Remarks/Alias	UNPH4; UNPH-2; USP15
Species	Human
GeneID (human)	9958
GeneID	9958
Immunogen	Synthetic peptide A synthetic peptide of human USP15
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4E8
Protein Weight	112 kDa
Shipping	Ice bag