


Anti-Usp19 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O94966

Item Number: CM0029240



Introduction

Rabbit polyclonal antibody specific for USP19 (ubiquitin carboxyl-terminal hydrolase 19), validated in WB, IF/ICC, ELISA. Cross-reacts with human and mouse; detects ~146 kDa band. Ideal for research on deubiquitination, ERAD, hypoxia, and muscle differentiation.

[Learn More](#)

Parameter	Value
Product Name	USP19 Rabbit pAb
Remarks/Alias	ZMYND9; USP19
Species	Human
GeneID (Human)	10869
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human USP19 (NP_001186089.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O94966
Protein Weight	146kDa
Shipping	Ice bag