


Anti-Usp21 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Uk80

Item Number: CM0024262



Introduction

Rabbit polyclonal antibody targeting human USP21 deubiquitinase, suitable for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. USP21 regulates histone H2A deubiquitination and transcriptional activation.

[Learn More](#)

Parameter	Value
Product Name	USP21 Rabbit pAb
Remarks/Alias	USP16; USP23; USP21
Species	Human
GeneID	27005
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human USP21 (NP_036607.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UK80
Protein Weight	63kDa
Shipping	Ice bag