

Anti-Usp22 Rabbit Monoclonal Antibody For Wb And Elisa - Q9Upt9

Item Number: CM0009255

Introduction

USP22 Rabbit Monoclonal Antibody validated for WB and ELISA, cross-reacting with human and mouse USP22. Ideal for studying deubiquitinase roles in transcription regulation, cell cycle progression, and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	USP22 Rabbit mAb
Remarks/Alias	USP3L; USP22
Species	Human
Gene ID	23326
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-175 of human USP22 (Q9UPT9)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UPT9
Protein Weight	60kDa
Shipping	Ice bag