

Anti-Usp22 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Upt9

Item Number: CM0018926

Introduction

USP22 Rabbit Polyclonal Antibody validated for WB and ELISA, reacting with human and rat.

USP22 is a deubiquitinase involved in transcriptional regulation, cell cycle, and innate immunity, suitable for research applications.

[Learn More](#)

Parameter	Value
Product Name	USP22 Rabbit pAb
Remarks/Alias	USP3L; USP22
Species	Human
GeneID (Human)	23326
GeneID	23326
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human USP22 (NP_056091.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9UPT9
Protein Weight	60kDa
Shipping	Ice bag