

Anti-Usp24 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Upu5

Item Number: CM0002355

Introduction

USP24 Rabbit monoclonal antibody validated for WB and ELISA. Detects human, mouse, and rat USP24. Ideal for ubiquitin-proteasome, cell survival, and ferritinophagy research.

[Learn More](#)

Parameter	Value
Product Name	USP24 Rabbit mAb
Remarks/Alias	USP24; ubiquitin specific peptidase 24
Species	Human
GeneID	23358
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 2500-2600 of human USP24 (Q9UPU5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UPU5
Protein Weight	294kDa
Shipping	Ice bag