

Anti-Usp24 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Upu5

Item Number: CM0020892

Introduction

Rabbit polyclonal antibody targeting human USP24, a 294kDa ubiquitin-specific protease involved in ferritinophagy and cell survival. Validated for Western blot and ELISA with human and mouse reactivity.

[Learn More](#)

Product Name	USP24 Rabbit pAb
Remarks/Alias	USP24; ubiquitin specific peptidase 24
Species	Human
Gene ID (Human)	23358
Gene ID	23358
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UPU5
Protein Weight	294kDa
Shipping	Ice bag