

Anti-Usp24 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Upu5

Item Number: CM0015784

Introduction

USP24 Rabbit Polyclonal Antibody (pAb) specifically detects human and mouse USP24, validated for Western blot and ELISA. This deubiquitinase antibody supports research on cell survival, ferritinophagy, and ubiquitin pathways.

[Learn More](#)

Parameter	Value
Product Name	USP24 Rabbit pAb
Remarks/Alias	USP24
Species	Human
GeneID (Human)	23358
GeneID	23358
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2431-2620 of human USP24 (NP_056121.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UPU5
Protein Weight	294kDa
Shipping	Ice bag