


Anti-Usp25 Rabbit Monoclonal Antibody For Wb - Q9Uhp3

Item Number: CM0005230



Introduction

Rabbit monoclonal antibody against human USP25, validated for WB, IP, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for deubiquitinase, inflammation, and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	USP25 Rabbit mAb
Remarks/Alias	USP21; USP25
Species	Human
GeneID	29761
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human USP25 (NP_037528.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UHP3
Protein Weight	122KDa/125KDa/129KDa
Shipping	Ice bag