


Anti-Usp25 (Usp21) Rabbit Polyclonal Antibody For Wb, Elisa - Q9Uhp3

Item Number: CM0017826



Introduction

High-quality rabbit polyclonal antibody against human USP25, a deubiquitinating enzyme involved in inflammation, Wnt signaling, and immune response. Validated for WB and ELISA, with cross-reactivity to mouse and rat.

[Learn More](#)

Parameter	Specification
Product Name	USP25 Rabbit pAb
Remarks/Alias	USP21; USP25
Species	Human
Gene ID (Human)	29761
Gene ID	29761
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human USP25 (NP_037528.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UHP3
Protein Weight	122kDa
Shipping	Ice bag