

Anti-Usp26 Polyclonal Antibody For Wb, Elisa - Q9Bxu7

Item Number: CM0024791

Introduction

USP26 Rabbit Polyclonal Antibody validated for WB and ELISA, reacts with human, mouse, and rat samples. Useful for studying deubiquitination, androgen receptor signaling, and spermatogenesis.

[Learn More](#)

Parameter	Value
Product Name	USP26 Rabbit pAb
Remarks/Alias	SPGFX6; USP26
Species	Human
Gene ID (Human)	83844
Gene ID	83844
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 664-913 of human USP26 (NP_114113.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BXU7
Protein Weight	104kDa
Shipping	Ice bag