

Anti-Usp30 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q70Cq3

Item Number: CM0013690



Introduction

Rabbit polyclonal antibody against human USP30, validated for Western blot, immunofluorescence, and ELISA. Reacts with human, mouse, and rat samples. Ideal for mitophagy and deubiquitination research.

[Learn More](#)

Parameter	Value
Product Name	USP30 Rabbit pAb
Remarks/Alias	USP30
Species	Human
GeneID (Human)	84749
GeneID	84749
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 57-390 of human USP30 (NP_116052.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q70CQ3
Protein Weight	59kDa
Shipping	Ice bag