


Anti-Usp33 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q8Tey7

Item Number: CM0019365



Introduction

High-quality rabbit polyclonal antibody targeting human USP33 (VDU1) for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Essential tool for deubiquitination, centrosome biology, and GPCR trafficking research.

[Learn More](#)

Parameter	Value
Product Name	USP33 Rabbit pAb
Remarks/Alias	VDU1; USP33
Species	Human
Gene ID	23032 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 290-400 of human USP33 (NP_963920.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TEY7
Protein Weight	107kDa
Shipping	Ice bag