

Anti-Usp39 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q53Gs9

Item Number: CM0003596



Introduction

Anti-USP39 rabbit monoclonal antibody validated for WB, IHC-P, ELISA. Reacts with human and rat. USP39 is a deubiquitinase involved in splicing, DNA repair, and cell cycle control. Ideal for research in B-cell development, DNA damage response, and antiviral signaling.

[Learn More](#)

Parameter	Value
Product Name	USP39 Rabbit mAb
Remarks/Alias	65K; SAD1; CGI-21; HSPC332; SNRNP65; USP39
Species	Human
GeneID (Human)	10713
GeneID	10713
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human USP39 (Q53GS9).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q53GS9
Protein Weight	65kDa
Shipping	Ice bag