


Anti-Usp39 Rabbit Polyclonal Antibody For Wb And Elisa - Q53Gs9

Item Number: CM0030796



Introduction

Rabbit anti-USP39 polyclonal antibody (alias 65K, SAD1, SNRNP65) validated for Western blot and ELISA. Shows reactivity with human, mouse, and rat USP39. Suitable for studying splicing, DNA repair, and antiviral signaling.

[Learn More](#)

Parameter	Specification
Product Name	USP39 Rabbit pAb
Remarks/Alias	65K; SAD1; CGI-21; HSPC332; SNRNP65; USP39
Species	Human
GeneID (human)	10713
GeneID	10713
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 336-565 of human USP39 (NP_006581.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q53GS9
Protein Weight	65kDa
Shipping	Ice bag