

## Anti-Usp4 Polyclonal Antibody For Wb, If/Icc, Elisa - Q13107

Item Number: CM0019679



### Introduction

Rabbit polyclonal antibody against human USP4, validated for WB, IF/ICC, and ELISA. Detects USP4 in human and mouse samples. Ideal for deubiquitination, immune response, and cancer biology studies.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | USP4 Rabbit pAb                                                                                                   |
| Remarks/Alias    | UNP; Unph; USP4                                                                                                   |
| Species          | Human                                                                                                             |
| GeneID (human)   | 7375                                                                                                              |
| GeneID           | 7375                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 560-700 of human USP4 (NP_003354.2) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, IF/ICC, ELISA                                                                                                 |
| Cross Reactivity | Human, Mouse                                                                                                      |
| SWISS            | Q13107                                                                                                            |
| Protein Weight   | 109kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                           |