

Anti-Usp47 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q96K76

Item Number: CM0016647

Introduction

USP47 rabbit polyclonal antibody validated for WB, IHC-P, ELISA detection of human and mouse USP47 (~157 kDa). Ideal for DNA repair, genome integrity, and cell cycle research. - Q96K76.

[Learn More](#)

Parameter	Value
Product Name	USP47 Rabbit pAb
Remarks/Alias	TRFP; USP47
Species	Human
Gene ID (Human)	55031
Gene ID	55031
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1128-1287 of human USP47 (NP_060414.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96K76
Protein Weight	157kDa
Shipping	Ice bag