

Anti-Usp48 Rabbit Polyclonal Antibody For Wb, Elisa - Q86Uv5

Item Number: CM0028713

Introduction

This USP48 rabbit polyclonal antibody recognizes human, mouse, and rat USP48 (deubiquitinase). Validated for WB and ELISA, it aids in studying NF-κB signaling, cell cycle, and pyroptosis regulation.

[Learn More](#)

Parameter	Value
Product Name	USP48 Rabbit pAb
Remarks/Alias	USP31, DFNA85, RAP1GA1, USP48
Species	Human
Gene ID (Human)	84196
Gene ID	84196
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human USP48 (NP_115612.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86UV5
Protein Weight	119kDa
Shipping	Ice bag