


Anti-Usp5 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P45974

Item Number: CM0026126



Introduction

USP5 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, rat USP5 (96 kDa). Ideal for deubiquitination and signaling pathway studies.

[Learn More](#)

Parameter	Value
Product Name	USP5 Rabbit pAb
Remarks/Alias	ISOT; USP5
Species	Human
GeneID (human)	8078
GeneID	8078
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human USP5 (NP_003472.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45974
Protein Weight	96kDa
Shipping	Ice bag