

Anti-Usp53 Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - Q70Ek8

Item Number: CM0021042

Introduction

Anti-USP53 Rabbit Polyclonal Antibody validated for WB, IF/ICC, IF-P, ELISA. Specifically detects human and mouse USP53, a 120 kDa tight junction deubiquitinase involved in hearing and barrier function.

[Learn More](#)

Parameter	Details
Product Name	USP53 Rabbit pAb
Remarks/Alias	USP53; ubiquitin specific peptidase 53
Species	Human
Gene ID (Human)	54532
Gene ID	54532
Immunogen	Recombinant protein corresponding to a sequence within amino acids 495-582 of human USP53 (NP_061923.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q70EK8
Protein Weight	120kDa
Shipping	Ice bag