


Anti-Usp8 Rabbit Pab For Wb, If/Icc, Elisa - P40818

Item Number: CM0029060



Introduction

USP8 Rabbit Polyclonal Antibody for Western blot, IF/ICC, and ELISA. Recognizes human, mouse, and rat USP8. USP8 is a deubiquitinase regulating protein turnover, endosomal sorting, and EGFR signaling. Ideal for cell cycle and proteostasis research.

[Learn More](#)

Parameter	Value
Product Name	USP8 Rabbit pAb
Remarks/Alias	UBPY; PITA4; SPG59; HumORF8; USP8
Species	Human
Gene ID	9101
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human USP8 (NP_005145.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40818
Protein Weight	128kDa
Shipping	Ice bag